



The Truth About Blood

ROM Plus is the ONLY immunoassay test that is quantifiably tested with a specified concentration of blood

10% blood concentration

ROM Plus

ROM Plus IFU - "Blood [was] spiked into the low positive at a 10% final concentration without loss of activity... [B]lood [was] also spiked into the negative-control matrix and shown to be negative."

"ROM Plus test kits will function properly with trace amounts of blood in the sample. Significant amounts of blood discharge may cause the test to malfunction and is not recommended."

Unknown blood concentration

Actim PROM™

Actim PROM IFU - "Whole blood with concentrations corresponding to typical pregnancy levels of IGFBP-1 was tested and did not affect Actim PROM test performance."

"The test has been designed to minimize interference from bleeding, but in cases of heavy bleeding the blood locally may have a higher concentration of IGFBP-1 protein. In these cases, a positive result should be interpreted with caution."

ActimPROM is a trademark of Medix Biochemica

Local Distributor:

clinicalinnovations.com.au
36 Matheson Street, Virginia
Queensland 4014 | Australia
Phone +07 3862 2000 | Fax +07 3862 2077
enquiries@clinicalinnovations.com.au



© 2020 Clinical Innovations, LLC. All rights reserved. | ART-0266 Rev. 01