



The Koala IU<sup>®</sup>PC is compatible with all major brands of fetal monitors. This tool is designed to help identify what fetal monitor is being used and to ensure the correct cables are ordered for an in-service or evaluation.

The most common cables that are used with the Koala are the **GE Medical Corometrics (CI IPC-5014) cable**, the **Philips (CI IPC-5080) cable** and the **Philips (CI IPC-5012) cable**.


CI IPC-5014

Compatible with **GE Corometrics 250cx**



CI IPC-5080

Compatible with **Philips Avalon Monitors**



FM50

FM30

As a reminder when zeroing the transducer while using the **Phillips Avalon monitor** the catheter needs to be disconnected from the cable, not the cable from the adapter cable.

CI IPC-5012

Compatible with  
**All Other Philips Monitors**





# KOALA CABLE TOOL

## CI IPC-5016

Compatible with **Spacelabs Medical AMS 94000**



## CI IPC-5090

Compatible with **EDAN Instruments**



The Koala Intrauterine Pressure Catheter is the only IUPC that has no electronics in the catheter.

The external transducer provides easy, true zeroing in utero.

### Ordering Information

Product Description	Product Code	Box Quantity	Monitor Cables	Product Code	Box Quantity
Koala Intrauterine Pressure Catheter	CI IPC-5000E	10/ea	Philips 8040A (40V/mmHg)	CI IPC-5018	1 Cable
			Philips Avalon	CI IPC-5080	1 Cable
			All other Philips Monitors	CI IPC-5012	1 Cable
			GE Medical (Corometrics) All Models	CI IPC-5014	1 Cable
			Air-Shields (Litton) Fetascan 1400	CI IPC-5016	1 Cable
			Spacelabs Medical (AMS) Spacelabs 94000	CI IPC-5050	1 Cable
			Neovent Medical Stan S31	CI IPC-5060	1 Cable
			Edan Instruments, Inc.	CI IPC-5090	1 Cable

See Instructions for Use for full instructions, warnings, precautions and contraindications

Local Distributor:

[clinicalinnovations.com.au](http://clinicalinnovations.com.au)

36 Matheson Street, Virginia  
Queensland 4014 | Australia  
Phone +07 3862 2000 | Fax +07 3862 2077  
enquiries@clinicalinnovations.com.au

